

Refrigerated Vaccine Inside at 2°C to 8°C



RECEIVING AND UNPACKING STEPS

- 1 IMMEDIATELY** open the cooler and follow the TagAlert receiver instructions. See below for details.
- 2 Compare** the packing list to the contents of the cooler to confirm they match. See below.
- 3 Inspect** the vaccine and report any tampering (e.g., broken seal). See back for details.
- 4 Place vaccine** in appropriate storage immediately upon receipt.
- 5 Repack** the cooler for refurbishing. See back for details.
- 6 Return** the cooler via UPS. See back for details.



INSTRUCTIONS FOR CONFIRMING PROPER SHIPPING TEMPERATURES

Receiver Instructions

- 1. Stop TagAlert® Indicator**
Press and hold the START & STOP button for 5 seconds until Stop icon ● appears.
- 2. Read Indicator Status**
✓ **Accept:** Left arrow ◀ points to Green check mark: Ready for Use.
✗ **Reject:** Right arrow ▶ points to Red X and numbers 1, 2 and/or 3 are visible—a temperature alarm has triggered. **QUARANTINE PRODUCT** and call +1-877-836-7123 or your state/local immunization program right away.

T17003135 Rev. A

TAGALERT INFORMATION

- Remove the TagAlert temperature monitor from the box. The TagAlert is located inside the inner box, next to the vaccine.
- Follow the recipient instructions on the blue card:
 - Stop TagAlert Indicator.** Press and hold the Start and Stop button until the Stop icon appears.
 - Read Indicator Status.** Determine if the appropriate temperature was maintained during transit.
 - TagAlert Monitor Concerns?** Please call **1-877-836-7123**.



Vaccine should be placed in appropriate storage immediately upon receipt.



REVIEW YOUR PACKING SLIP

McKesson Specialty Care Dist.
8745 Forest Hill Irene, Suite 105
Olive Branch, MS 38654

Packing Slip

This is not an invoice

Page 1 of 2

Ship-to:

McKesson Specialty Care Dist.
8745 Forest Hill Irene, Suite 105
Olive Branch, MS 38654

Awarde

McKesson Specialty Care Dist.
8745 Forest Hill Irene, Suite 105
Olive Branch, MS 38654

Provider PIN: WVA200066

Delivery Number: 230684533

Quality Check Date: 02/11/2025

Internal use only



NDC	Customer P.O.	Material Description Manufacturer	MFR Lot#	Exp. Date	VFC Doses	317 Doses	State Doses	CHIP Doses	Order Qty	Ship Qty	Unit Price	Extended Price
19515-0810-52 5017406	394261	FLU; SYR; 10-pack (TIV) PED CDC GLAXO SMITHKLINE	AE2J7	06/16/2025	10	0	0	0	10	10	\$15.22	\$152.20
Total									10	10		\$152.20

Verify the **MFR Lot#** and the **Ship Qty** on the packing slip with the contents of the cooler.
The **Ship Qty** is the total number of vaccine cartons.

Refrigerated Vaccine Inside at 2°C to 8°C



RETURN YOUR ECOFLEX COOLER

**ECOFLEX
MUST BE
RETURNED**
Do Not Dispose



Step 1: Open shipper and remove the top PCM gel pack(s).



Step 2: Remove the product(s) from the inner corrugated box.



Step 3: Ensure the inner corrugated box is back in the EcoFlex shipping system.



Step 4: Place the top PCM gel pack(s) back on top.



Step 5: Fold in the original outer flaps. Fold in the remaining flaps with return shipping label exposed. Tape the seam.



Step 6: UPS will arrive within 1 – 2 days for EcoFlex pickup.



UPS WILL AUTOMATICALLY PICKUP ECOFLEX COOLERS



- 1 UPS will arrive within 1 – 2 days after delivery for EcoFlex pickup.
- 2 When the driver arrives, return all EcoFlex coolers.



DO NOT RETURN EPS COOLERS



**Do NOT Return KoolTemp Coolers
PLEASE DISPOSE**

**DO NOT
RETURN
KOOLTEMP
COOLERS**
Please Dispose